



Product Technical Information

Effective: May 27, 2003
Re: Latex Content

The Agee Digit Widget™ External Fixation System is an external fixator for treatment of proximal interphalangeal joint flexion contractures.

There are elastic bands in the device which contain natural rubber latex which may cause allergic reactions.

These 1/4" diameter elastic bands are not part of the surgical portion of the kit. They are applied by medical staff as part of the non-surgical assembly of the device then changed by the patient on a daily basis. These bands are similar to bands used for orthodontic braces.

There are no other latex components in the Digit Widget. Contact Hand Biomechanics Lab, Inc. at 800-522-5778 if you need further technical information on the Digit Widget System.