

**\*\* BULLETIN \*\***

**Hand Biomechanics Lab starts  
REMANUFACTURING program  
for our AGEE WRISTJACK® System**


Hand Biomechanics Lab, Inc., the original manufacturer of the Agee WristJack, has just received FDA's approval of our 510(k) application to market remanufactured WristJack Systems. Participating in this program is an easy way to reduce equipment costs while still maintaining the high quality you've come to expect from us. We remanufacture your used WristJacks and return them directly to you.

To participate in our program, call and request a return kit. This kit contains detailed instructions, packaging and postage needed to return both the instrumentation and device. The instrumentation is conveniently returned right after the case. When the device is removed from the patient (usually 6-8 weeks later), obtain it from the surgeon, call us with a purchase order, then send us the WristJack according to the return kit instructions.

Once received by Hand Biomechanics Lab, we will disassemble, clean and inspect the returned WristJack and instruments. Your device is reassembled, replacing wear sensitive parts as needed, then fully tested to the same standards as a new device. The WristJack, instruments, new fixation pins and new surgeon's manual are repackaged in NON-STERILE form and sent back to you. You will need to steam sterilize them before use. Allow about 4 weeks for return of your remanufactured WristJack System, which remains a SINGLE USE device.

Catalog No: *CFD-147-RNS*  
Product Name: *Agee WristJack External Fixation System, REMANUFACTURED*  
Price: *\$825*

Note: The maximum number of times a new WristJack can be remanufactured is twice, giving you a total of three patient uses. Look for an embossed indicator indicating the number of times a WristJack has been remanufactured.



Customer Service: 800-522-5778

Remanufacturing Information Hotline: 888-974-7852